

VISCOUNT POOLS of Canton

PRO / VISION PURIFICATION SYSTEM

(PREMIUM LOW CHLORINE ALTERNATIVE)

New Season Start Up

1. After the pool is filled and the filter system is operating properly. Check and adjust the water chemistry (PH, Alkalinity, Calcium and Phosphates) to ensure the new chemicals work correctly, maintain quality swimmer comfort level and to safeguard equipment longevity and warranties remain valid.

(Free water testing at Viscount of Canton)

(PH: 7.2 - 7.8, Alkalinity: 80 - 120, Calcium Hardness: 200 - 300, Phosphates: 0 - 500)

2. Test and adjust the waters Cyanuric Acid (**Stabilizer or Conditioner**) to stabilize Chlorine.

(Cyanuric Acid: 30 - 75) * Slow dissolving chemical, put in a nylon or sock and hang in front of the pools return jet. This may take a few weeks to dissolve in cold water.

3. Add _____ lbs. **Pool Shock** (1lb. per 10'000 gallons of fresh water) Circulate 24 hours
(If water is cloudy or discolored Double Shock with 2lbs. per 10'000 gallons)

4. Add _____ lbs. **Pro Team Supreme Plus** (20lbs. Per 10'000 gallons)

***** Pro Team Supreme Plus does not leave the pool unless water is splashed out or drained, test level before adding starting amount each new season *****

(Install a new VISION Nature 2/9 Limited CRT Mineral Cartridge to start each new season)

Weekly Maintenance

Wednesday:

1. Check Water Chemistry and adjust if necessary .
2. Add _____ cups **Viscount Premium Soft Cal Chlor Ultra AC Granular** (1 cup per 5'000 Gallons)

***** Viscount Quick Tabs should ONLY be placed in the skimmer basket with filter on high speed until dissolved. Never put Quick Tabs in a Floater or Chlorinator *****

Sunday:

1. Add _____ cups **Viscount Premium Soft Cal Chlor Ultra AC granular** (1 cup per 5'000 Gallons)
2. Add _____ oz. **Viscount Nature Clear Enzyme Treatment** (4oz per 10'000 Gallons)

***** Viscount of Canton recommends to Filter / Circulate your pool water on High Speed 12 hours daily (night time filtration is more effective and cost efficient) or 24 hours on Low Speed *****